HYPERHYDROPHILIC

hydrorise[®]

VINYLPOLISILOXANE (ADDITION SILICONE)
IMPRESSION MATERIAL

Putty

Heavy Body

Monophase

Regular Body

Light Body









Extra Light Body

Zhermack[®]
CLINICAL 





Zhermack SpA
Via Bovazecchino, 100
45021 Badia Polesine (Rovigo) - Italy
tel. +39 0425 597611 - fax +39 0425 53596
www.zhermack.com - info@zhermack.com

Hydrorise Maxi 5:1 (380 ml)
Instruction for use with Modulmix,
Zhermack Automatic Mixing machine

Hydrorise 1:1 (50 ml)
Instruction for use with D2 dispenser

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FI EL DA SV SL RU PL

CS TR RO HR HU SK LV

LT ET CN AR BG JP KO

Table 1. Technical Data

	Hydrorise	Putty		Heavy Body	
1	ISO 4823	Type 0 Putty Consistency		Type 1 Heavy-bodied Consistency	
2	ADA Spec. 19	Type 0 Putty Consistency		Type 1 Heavy-bodied Consistency	
3	Mixing ratio (Base:Catalyst)	1:1 5:1		1:1 5:1	
4	Strain in Compression (Min.–Max.)	1 – 3%		3 – 5%	
5	Recovery from deformation	>99,0%		>99,5%	
6	Linear dimensional change (24hrs)	<0,20%		<0,20%	
7	Hardness "Shore A"	60±2		60±2	
	TIMES FOR CLINICAL USE*	Normal Setting	Fast Setting	Normal Setting	Fast Setting
8	Mixing Time	30" (1:1) Auto mixing System		Auto mixing System	
9	Total Working Time	2'00"	1'30"	2'00"	1'30"
10	Time in Mouth (minimum)	3'30"	2'30"	3'30"	2'30"
11	Total Setting Time	5'30"	4'00"	5'30"	4'00"
12	Aroma and Flavour	The products in the Hydrorise range have a vanilla mint aroma			

*The times mentioned must be intended from the start of mixing at 23°C – 73°F.

Monophase		Regular Body		Light Body		Extra Light Body	
Type 2 Medium-bodied Consistency		Type 2 Medium-bodied Consistency		Type 3 Light-bodied Consistency		Type 3 Light-bodied Consistency	
Type 2 Medium-bodied Consistency		Type 2 Medium-bodied Consistency		Type 3 Light-bodied Consistency		Type 3 Light-bodied Consistency	
1:1 5:1		1:1		1:1		1:1	
3 – 5%		3 – 5%		3 – 5%		3 – 5%	
>99,5%		>99,5%		>99,5%		>99,5%	
<0,20%		<0,20%		<0,20%		<0,20%	
54±2		45±2		45±2		45±2	
Normal Setting	Fast Setting	Normal Setting	Fast Setting	Normal Setting	Fast Setting	Normal Setting	Fast Setting
Auto mixing System		Auto mixing System		Auto mixing System		Auto mixing System	
2'00"	1'30"	2'00"	1'30"	2'00"	1'30"	2'00"	1'30"
3'30"	2'30"	3'30"	2'30"	3'30"	2'30"	3'30"	2'30"
5'30"	4'00"	5'30"	4'00"	5'30"	4'00"	5'30"	4'00"
and flavour							

Higher temperatures reduce the times, lower temperatures increase them.

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SOLO PER USO DENTALE

HYDROPHILIC ADDITION SILICONES (VINYL POLYSILOXANE) FOR HIGH-PRECISION IMPRESSIONS

USE

Hydrorise Putty: Vinyl polysiloxane for first impressions with elastic final consistency, recommended for all impression techniques: two-stage(double impression*) or single-stage (simultaneous).

Hydrorise Heavy Body: High viscosity vinyl polysiloxane. Recommended for filling impression trays in conjunction with a 'wash' within the Hydrorise range (i.e. light body, extra light body & regular body) in the single-stage (simultaneous) technique.

Hydrorise Monophase: Medium viscosity hyper-hydrophilic vinyl polysiloxane. A single viscosity monophase material recommended for use as a preparation wash and for filling individual impression trays. Also if desired this Monophase can be used in conjunction with a 'wash' from the HYdrorise range (i.e. light body, extra light body & regular body) in the single-stage (simultaneous) technique.

Hydrorise Regular Body: Medium viscosity hyper-hydrophilic vinyl polysiloxane. Recommended for use in conjunction with Hydrorise Putty or Heavy Body in the single-stage (simultaneous) technique, or as a monophase material in a special tray.

Hydrorise Light Body: Low viscosity hyper-hydrophilic vinyl polysiloxane. Recommended for use in conjunction with Hydrorise Putty or Heavy Body in the two-stage (double impression*) or single-stage (simultaneous) technique.

Hydrorise Extra Light Body: Very low viscosity hyper-hydrophilic vinyl polysiloxane. Recommended for use in conjunction with Hydrorise Putty or Heavy Body in the two-stage (double impression*) or in the single-stage (simultaneous) technique.

*When using this product in the two-stage (double impression) technique, we recommend that you place a polythene spacer sheet on the putty during the first stage of impression taking to improve detail quality in the next stage.

IMPRESSION TRAYS: PREPARATION AND ADHESIVES

All standard and/or special trays can be used, provided they are rigid and non-deformable. The range of **Zhermack's Hi-Tray** impression trays are especially recommended.

Proper adhesion of the material to the impression tray is essential to avoid errors, especially during removal from the mouth. To ensure a strong bond, make sure that the adhesive is appropriate for the type and material of the tray. Three different types of adhesive are supplied for this purpose:

IPERLINK LCT: A two-component adhesive for light-curing special trays.

ELITE IPERLINK SCT: An adhesive for self-curing special trays and disposable plastic trays.

UNIVERSAL TRAY ADHESIVE: An especially effective adhesive for standard metal trays.

IMPORTANT:

When using special trays, always ask the dental technician whether the tray is self-curing or light-curing.

For correct use, follow the instructions for each individual adhesive.

Hydrorise Putty (mixing ratio 1:1)

MIXING

- Take the required amounts of both components out of the tubs using the colour coded dosing scoops.

- Mix identical proportions of Hydorise Putty Base and Catalyst until you obtain a homogeneous mix with no streaks in it (mixing time: 30 seconds).

WARNINGS

- Respect the specified proportions precisely to ensure correct setting times. A higher dosage of catalyst does not speed up setting.
- Close the tubs immediately after use. Take care not to mix up the lids or dosing scoops.
- Contact with latex gloves may inhibit vinyl polysiloxane from setting. Avoid direct or indirect contact with latex gloves. Wash your hands thoroughly to eliminate all traces of impurities or use of vinyl or nitrile gloves.
- Vinyl polysiloxane is chemically resistant, so take care not to stain your clothing.

Hydorise Putty, Heavy Body, Monophase (Maxi Version 380ml, 5:1)

PREPARING THE DEVICE

- Press in the levers and remove the protective cap (fig. 1).
- Insert the cartridge into the machine as instructed by the manufacturer.
- Ensure an equal length of material is coming out. Clean and remove any excess material (fig. 2).
- Fit and secure the dynamic-static mixing tip (fig. 3).
- Fit the **mixing tip locker** taking care to engage the clips (fig. 4A) on the cartridge (fig. 4B).
- The Zhermack 5:1 system cartridge is now ready to use (fig. 6).
- Leave the dynamic-static mixing tip in place after use.
To replace it, simply turn the **mixing tip locker** through 45° and remove it (fig. 7A). - Then remove the dynamic-static mixing tip (fig. 7B).
- Follow the machine manufacturer's instructions to remove and replace the cartridge.

WARNINGS

- Always dispense equal amounts of base & catalyst before fitting the dynamic-static mixing tip. If you forget to do so, start the machine and wait until the material coming out is uniform in colour, then take away any non-homogenous material before proceeding.
- Only use original Zhermack 5:1 system cartridges, dynamic-static mixing tips and **mixing tip lockers**. The use of other mixing tips does not guarantee optimal mixing and can jeopardise the product's characteristics and end results.
- The **mixing tip locker** cannot be used with cartridges not from Zhermack's 5:1 mixing system.

Hydorise Light Body, Extra Light Body, Regular Body (cartridge version 50ml, 1:1)

PREPARING THE DEVICE

- Assemble the device as illustrated in the diagrams (figs. 8-9).
- Before fitting the tip, make sure that the two components (base and catalyst) flow out evenly by pressing gently on the dispenser lever and extruding a small amount of material. Remove this material before proceeding (fig. 10).
- Fit the mixing tip and optional intra-oral tip on to the cartridge (fig. 11).

MATERIAL CHARACTERISTICS - (Technical Data)

Table 1 gives the characteristics of the materials:

Row 1 and 2: Consistency: ISO 4823 and ADA Spec 19 respectively.

Row 3: Mixing ratio for base and catalyst

Row 4: Strain in compression

Row 5: Recovery from deformation

Row 6: Linear dimensional change at 24 hours

Row 7: Shore-A hardness

Times for clinical use*

Row 8: Mixing time

Row 9: Total working time

Row 10: Time in mouth

Row 11: Total setting time

Row 12: The products in the Hydrorise range have a vanilla mint aroma and flavour: see Table 1

*The times for clinical use (listed in **table 1 – Table 1**, material characteristics) are intended from the start of mixing at a temperature of 23°C – 73°F. Mixing at higher temperatures reduces these times. Mixing at lower temperatures increases them.

WARNING:

In sensitive subjects, vinyl polysiloxane may cause irritation or other allergic reactions.

CLEANING AND DISINFECTING THE IMPRESSION

After rinsing well under running water, the impression can be immediately disinfected. Ideal disinfecting is achieved by immersing the impression in **Zeta 7 solution** or spraying it with **Zeta 7 spray** (Zhermack). When using other disinfectants, always follow the manufacturer's instructions.

KEEPING IMPRESSIONS AND CASTING MODELS

There is no need to wait. Models may be cast immediately after disinfection or up to two weeks after, provided the impression is kept at room temperature. Hydrorise is compatible with all the best stones/plasters on the market and with polyurethane resins. We recommend **Elite Rock / Elite Rock Fast** (Zhermack) extra hard die stone (class 4) and **Elite Model / Elite Model Fast** (Zhermack) for hard plaster models (class 3).

CLEANING AND DISINFECTING THE DISPENSER

To clean, disinfect or sterilize the dispenser, always follow the manufacturer's instructions. We recommend **Zeta 3** (Zhermack) disinfectant and **Zeta 2** (Zhermack) sterilisation liquid.

STORAGE AND WARRANTY

Hydrorise impression materials are guaranteed for 36 months from the date of production provided they are stored correctly between 5°- 27°C / 41°- 80°F.

GALVANISATION

Hydrorise impressions can be galvanised with copper or silver.

IMPORTANT REMARKS: Consulting services on the use of our products, provided in any form, be that verbal, written or through demonstrations, are based upon the current state of dentistry and of our know-how. They are to be considered as non-binding information, also in relation to possible third party claims, and do not exonerate the user from personally checking whether the product is suitable for the planned use. Use and application by the user take place under no possible control by the Company, and are therefore under the user's responsibility. Possible liability for damages is limited to the value of goods supplied by the Company and utilized by the user.

FOR DENTAL USE ONLY