

discovery[®]
smart



Small and brilliant.

D
DENTAURUM

The finishing touches make the difference!



**Significant
reduction in size.**

20 % smaller than
discovery® and one of
the smallest brackets
in the world.



For over 20 years discovery® has stood for perfect design and maximum quality in bracket treatment. This proven premium bracket family now has a new member named discovery® smart. A bracket that is esthetically small but impressively large in features & function.

A practical feature of discovery® smart is its mesial-distal curved slot contour, which follows the ideal contour of the dental arch. This means that each bracket has an individual slot curvature, which helps to minimize the friction in the slot and ensures more efficient & precise treatment.

The dimensions of this new family member are 20 % smaller than discovery®.



Outstanding sliding properties.

New mesial-distal curved slot contour aligned to the ideal arch.

Brilliant properties.



One-piece bracket.

High biocompatibility.

Optimal wing design.

Easy, efficient ligating.

The bracket is manufactured as a one-piece bracket on state-of-the-art equipment using the MIM technique (metal injection molding), and has a very low profile with an in-out value of only 0.7 mm. The perfectly designed tie wings, provide the bracket with the perfect finish.

The small, modern design of discovery® smart also incorporates all the brilliant characteristics of the discovery® bracket range, ensuring maximum reliability and control with fascinating esthetics.




Break point.

Flex for reliable debonding.

Laser-structured base.

Patented laser-structured base for better bond and FDI notation for ease of use.



Mushroom-design hooks.

Secure retention of attachments.

Laser marking.

Clear, easy identification.

Low profile.

Excellent patient comfort.

3D bracket base.

Perfect fit on the tooth.

Multifaceted treatment.



In addition to discovery® smart, the discovery® family offers several other bracket choices which fulfill all requirements of a variety of treatment protocols and delivers a diverse combination of options within the range.

Brackets for all types of treatment.

discovery[®]
smart

Small and brilliant.





discovery[®]

The sum of experience.



discovery[®] *sl* *discovery*[®]
sl 2.0

High-tech for your success.





www.dentaurum.de

➔ For more information on our products and services, please visit www.dentaurum.de

Date of information: 06/12

Subject to modifications



D
DENTAURUM

Turnstraße 31 | 75228 Ispringen | Germany | Phone +49 72 31/803-0 | Fax +49 72 31/803-295
www.dentaurum.de | info@dentaurum.de